

Master Lot Sheet

Description

Product Name	TruSeq® Custom Amplicon Kit Dx		
Part Number	20005718	Catalog Number	20005718
Lot Number	A135913	Kit Expiration Date	10-JUL-2020

Master Kit Contents

Box 1 of 3

Item Description	Part Number	Lot Number
TRUSEQ® CUSTOM AMPLICON KIT DX 1/3	20003227	A135913-1
TRUSEQ® CUSTOM AMPLICON DX 1A	20003228	A135913-11
DX#-OHS1, Hybridization Buffer	15035442	20408510
DX#-TDP1, PCR Polymerase	15035469	20379640
DX#-PMM2, PCR Master Mix	15035448	20393893
DX#-A501, Index Primer A	15035449	20379611
DX#-A502, Index Primer B	15035450	20379612
DX#-A503, Index Primer C	15035451	20379613
DX#-A504, Index Primer D	15035452	20379648
DX#-A505, Index Primer E	15035453	20379614
DX#-A506, Index Primer F	15035454	20344837
DX#-A507, Index Primer G	15035455	20393806
DX#-A508, Index Primer H	15035456	20393807
DX#-A701, Index Primer 1	15035457	20379619
DX#-A702, Index Primer 2	15035458	20393809
DX#-A703, Index Primer 3	15035459	20379622
DX#-A704, Index Primer 4	15035460	20393811
DX#-A705, Index Primer 5	15035461	20393812
DX#-A706, Index Primer 6	15035462	20393813
DX#-A707, Index Primer 7	15035463	20403152
DX#-A708, Index Primer 8	15035464	20379629
DX#-A709, Index Primer 9	15035465	20403156
DX#-A710, Index Primer 10	15035466	20379645
DX#-A711, Index Primer 11	15035467	20379631

Item Description	Part Number	Lot Number
DX#-A712, Index Primer 12	15035468	20379634
DX#-ELM3, Extension – Ligation Mix	20005678	20423124
TRUSEQ® CUSTOM AMPLICON DX 1B	20003229	A135913-12
DX#-LNA1, Library Normalization Diluent	15035473	20379646
DX#-PX3, PhiX Internal Control	15036574	20393881
DX#-HT1, Library Dilution Buffer	20004095	20402886

Box 2 of 3

Item Description	Part Number	Lot Number
TRUSEQ® CUSTOM AMPLICON KIT DX 2/3	20003460	A135913-2
DX#-LNS1, Library Storage Buffer	20004094	20362530
DX#-EBT, Elution Buffer	20004093	20390433
DX#-FPU, Filter Plate	20012590	20389017

Box 3 of 3

Item Description	Part Number	Lot Number
TRUSEQ® CUSTOM AMPLICON KIT DX 3/3	20003463	A135913-3
TRUSEQ® CUSTOM AMPLICON DX 3A	20003461	A135913-31
DX#-SW1, Stringent Wash Buffer	15035443	20408513
DX#-UB1, Universal Wash Buffer	20004092	20388005
TRUSEQ® CUSTOM AMPLICON DX 3B	20003462	A135913-32
DX#-XPB, PCR Clean Up Beads	15036674	20409700
DX#-LNW1, Library Normalization Wash	15035446	20422776
DX#-LNB1, Library Beads	20017098	20379647